



S.R.L SUB – PRODUCTOS GANADEROS

República del Líbano 4279 - Quilmes - Bs. As. Argentina TEL/Fax: (011) 3551-5469

[infocanadas@hotmail.com](mailto:infocanadas@hotmail.com)

---

## **MEMORIA DESCRIPTIVA DEL PROCESO DE OREJAS DE CERDO DESHIDRATADAS PARA MASCOTAS**

- 1.- Se recibe la materia prima de establecimientos elaboradores registrados en SENASA.
- 2.- Se verifica la calidad primeramente estableciendo un parámetro de las calidades, con un análisis organoléptico de la materia prima.
- 3.- Se cargan las orejas a fulones procediéndose al lavado primario con agua para eliminar la suciedad que pueda venir en las mismas.
- 4.- Se escurre y se agregan salmonecidas para asegurar la sanitización del producto.
- 5.-Se recorta para eliminar restos de grasa que impiden el secado deseado.
- 6.- El producto se extiende en bandejas que a su vez se colocan en carros y es llevado a los hornos de deshidratación.
- 7.- Sobre el producto seco se realiza una primera clasificación de calidad y luego a una segunda clasificación de tamaño.
- 8.- Se realiza una clasificación definitiva y se sacan muestras para laboratorio (este paso es tercerizado a un laboratorio aprobado por SENASA)
- 9.-Ya clasificadas, las orejas se cuentan y se envasan de acuerdo al requerimiento del cliente.

## **DEHYDRATED PIG EARS DESCRIPTIVE MEMORY PROCESS**

- 1) Raw material is received from SENASA certificated establishments.
- 2) Raw material quality is verified based on a specific parameter using an organoleptic analysis.
- 3) The ears go into a wood fulon to proceed with the primary cleaning to eliminate the surface dirt with water.
- 4) The product is drained and sanitized with salmonelicide.
- 5) The ears edges are cut to eliminate left over fat allowing a desire drying.
- 6) The product is spread on trays that are placed on carts and taken to dehydration ovens.
- 7) A first quality classification is made on the dry product, and then a second size classification.
- 8) A definitive classification is made and samples are taken to a laboratory (this step is outsourced to a SENASA certified laboratory).
- 9) Once classified, the ears are counted and packed according to the customer's requirement

### **Presentación**

BOLSA POR 300 UNIDADES

PESO NETO: 9 KG

PESO BRUTO: 9,1 KG